

## Liko Silhouette Sling (Mod. 22)



- Pliable sling that adapts to the body and requires minimal wheelchair space
- Recommended for lifting to sitting to/from close-fitting wheelchairs
- Provides a slightly reclined sitting posture suitable for patients with poor trunk stability
- Adjustable head support
- Suitable to be used during bathing
- Available in Mesh

### Warranty Information

- 1 Year



Contract approved

**ndis**

**Item Code:**

05\_123603164\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

**Compliance**

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394450	Silhouette Mod 22 Sling - Mesh - XSMALL	200 kg	Mesh (Loops)
LSS394460	Silhouette Mod 22 Sling - Mesh - SMALL	-	-
LSS394470	Silhouette Mod 22 Sling - Mesh - MEDIUM	-	-
LSS394440	Silhouette Mod 22 Sling - Mesh - LARGE	-	-