

Supergrip Cutlery



- The Supergrip Cutlery is a comfortable, wide, built-up ribbed handled utensil set ideal for users with limited hand function
- Stainless steel, dishwasher safe
- Handle Length: 102 mm
- Handle Diameter: 38 mm
- Weight: 45 - 55 g

Warranty Information

- 1 Year

- ✔ Contracted DVA product
- ✔ Contract approved

Support @ Home

Function: Eating and drinking

Category: Low risk

Click [here](#) for more information

NDIS

Item Code:

05_150900319_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Unit Weight
DLK275310	Teaspoon	45 g
DLK275320	Tablespoon	55 g
DLK275330	Fork	45 g

Code	Display Name	Unit Weight
DLK275340	Knife	55 g