

Supergrip Cutlery



• The Supergrip Cutlery is a comfortable, wide, built-up ribbed handled utensil set ideal for users with limited hand function

• Stainless steel, dishwasher safe

Handle Length: 102 mmHandle Diameter: 38 mm

• Weight: 45 - 55 g

Warranty Information

• 1 Year	
----------	--

\bigcirc	Con	tracte	ed D'	VA p	rodu	ıct

Ocontract approved

ndis

Item Code:

05_150900319_0103_1_2 Click <u>here</u> for more information

Compliance

No applicable compliance standard

Code	Display Name
DLK275310	Teaspoon
DLK275320	Tablespoon
DLK275330	Fork
DLK275340	Knife