

## DermaSaver Leg Tube Protection



- Protects fragile skin from pressure, friction and minor traumas that result in tears, bruising and minor abrasion
- The Dermasaver garments are soft and loose with ample air-flow reaching the skin constructed of two layers of soft, lightweight microspring textile, which allows air-flow and moisture transfer
- Anti-microbial - will not support growth of bacteria, latex free totally washable
- Sizing: Simply measure the circumference at the widest part of the calf. For borderline measurements, choose the larger size for extra comfort

- ✓ Contracted DVA product
- ✓ Contract approved

Code	Display Name	Calf/Arm Circumference
PTA533317	DermaSaver - Leg Tube - X Sml	18-23cm
PTA533317L	DermaSaver - Leg Tube - Lge	38-48cm
PTA533317M	DermaSaver - Leg Tube - Med	30-38cm
PTA533317S	DermaSaver - Leg Tube - Sml	23-30cm
PTA533317XL	DermaSaver - Leg Tube - XL	48-58cm
PTA533318	DermaSaver - Leg Tube with Double Knee - X Sml	18-23cm
PTA533318L	DermaSaver - Leg Tube with Double Knee - Lge	38-48cm

Code	Display Name	Calf/Arm Circumference
PTA533318M	DermaSaver - Leg Tube with Double Knee - Med	30-38cm
PTA533318S	DermaSaver - Leg Tube with Double Knee - Sml	23-30cm
PTA533318XL	DermaSaver - Leg Tube with Double Knee - XL	48-58cm