

## Meiko Pan & Bottle Rack



- Drying and storage rack for care utensils and urine bottles
- Stainless steel
- Rack at the bottom
- Drip tray is sold separately, size is compatible with storage rack

 Contract approved

### Support @ Home

**Function:** Collecting urine and faeces

**Category:** Prescribed

Click [here](#) for more information



Code	Display Name	Total Length	Total Width	Total Height	Capacity
UPA400250	Meiko Pan & Bottle Rack - TUS140 - 9 Pans & 9 Bottles	350 mm	1400 mm	650 mm	9 Pans & 9 Bottles
UPA400260	Meiko Pan & Bottle Rack - TUS50 - 6 Pans & 6 Bottles	350 mm	500 mm	650 mm	6 Pans & 6 Bottles