

Gel Armrest Trough Liner

- Fits onto the Ottobock trough channel forearm pad armrest
- Provides a comfortable, supportive gel support pad for the user



Warranty Information

- 1 Year

✔ Contract approved

Support @ Home

Function: Body-worn assistive products for body protection

Category: Under advice

Click [here](#) for more information

Compliance

No applicable compliance standard



Code	Display Name	Total Length
PTA508470	Gel Armrest Trough Liner - Medium	330.2 mm
PTA508480	Gel Armrest Trough Liner - Large	482.6 mm