

Aspire Toilet Seat Raiser with Lid



- Moulded plastic seat raiser with plastic lid to mask odours
- Fits all standard toilet seats
- Adjustable mounting brackets to allow for correct fitting to different size toilets
- Available in three different height versions to suit individual requirements
- SWL: 200 kg

Warranty Information

- 1 Year



✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_091218063_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Material	Unit Weight
BTT149900	Aspire Toilet Seat Raiser with Lid- 50mm	Plastic	1.4 kg
BTT150100	Aspire Toilet Seat Raiser with Lid - 100mm	Plastic	1.4 kg

Code	Display Name	Material	Unit Weight
BTT150400	Aspire Toilet Seat Raiser with Lid - 150mm	Plastic	1.4 kg