

Coway BA-13 Bidet



The Coway BA13 Bidet is a luxurious, sleek, and versatile bidet designed to cater to the entire family. With Coway's long-standing expertise in toilet hygiene and water purification, this bidet is the pinnacle of personal health and hygiene solutions. From its self-sterilising, stainless steel twin spray nozzles to its heated dryer, the BA13 is packed with high-tech features for a superior cleansing and comfort experience.

- **Remote Control:** Easy-to-use, handheld or wall-mounted control for all functions.
- **Heated Features:** Heated dryer, seat, and water with adjustable temperature settings.
- **Massage/Move:** Rhythmic nozzle movement for thorough cleansing.
- **Feminine and Rear Cleansing:** Separate functions and nozzles for enhanced hygiene.
- **Self-Cleaning Nozzles:** Automatic nozzle cleaning before and after each use.
- **Quick Release:** Easily removable for cleaning and fitting to various toilet types.

The handheld remote reduces the need for reaching, minimising injury risk. It adjusts nozzle position, temperature, and water pressure.

Warranty Information

- 2 Years Limited - 1 Year Standard Manufacturer's Warranty, 1 Year Complimentary Parts & Labour Warranty

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_093318082_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Total Length	Total Width	Total Height	Unit Weight
BTT150120	Coway BA-13 Bidet - Warm Water - Elongated	528 mm	382 mm	144 mm	6 kg - 6.5 kg
BTT150140	Standard Electric Bidet	493 mm	382 mm	144 mm	6 kg - 6.5 kg