

ROHO Select Series Cushion Cover



permobil

ROHO

- This is available in Heavy Duty Cushion Cover, a fluid-resistant, textile enclosure.
- Replacement Covers for the following Cushions:
 - QUADTRO SELECT[®] HIGH PROFILE[®] Cushion
 - QUADTRO SELECT[®] MID PROFILE[®] Cushion

 Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name
PTC597800	ROHO Select Series Cushion Cover - 16.75x16.75in
PTC597810	ROHO Select Series Cushion Cover - 18x18in
PTC597910	ROHO Select Series Cushion Cover - 22x20in