

Gel Knee Protector Buttons

- Provide a soft surface to protect the users knees from bumping against legrests
- Attach to the inside of the legrests on a wheelchair
- Includes mounting hardware



Warranty Information

- 1 Year

 Contract approved

Support @ Home

Function: Accessories for wheelchair seating

Category: Prescribed

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name
PTA508800	Gel Knee Protector Buttons - SA Legrest - 3.18cm
PTA508820	Gel Knee Protector Buttons - Elevating Legrest - 3.18cm
PTA508810	Gel Knee Protector Buttons - SA Legrest - 1.98cm
PTA508830	Gel Knee Protector Buttons - Elevating Legrest - 1.98cm