

Aspire Raphael X Custom and X/X5 Classic Backrest



Customise the Aspire Raphael X Lift Chair to match your unique style and house
Available in a range of commercial quality fabric for long lasting durability

- **Soft Touch Fabric** - The soft touch weave and diverse colour palette makes the Manisa range a highly practical and functional solution These fabrics are comforting, easy to clean and durable, yet still soft to touch
- **Leather** - A luxurious artisan leather with a soft glove-like texture, will retain its unique suppleness, provided they are regularly cared for. The leather is produced and tested to international standards for properties such as rubbing, flex resistance, light fastness, thickness, and strength
- **Microsuede** - A Suede like texture combined with high performance characteristics ensures the Eastwood fabric range remains easy to clean and maintain
- **Healthcare Vinyl** The Lustrell Ranger range is high performance vinyl with HealthGuard Antibacterial and Antimicrobial treatment to assist with infection control and hygiene
- Custom order lead time varies, approximately 5 - 6 weeks depending on fabric availability
- Custom orders are subject to terms and conditions, please see page two on the order form below
- All colour swatches shown may not be a true representation of the actual product and should be used as a guide only
- All colour swatches shown may not be a true representation of the actual product and should be used as a guide only
- The Classic Backrest is available as an additional accessory for users who prefer a more traditional support profile for the Raphael X or X5 Lift Recline Chair

 Contract approved

Support @ Home

Function: Accessories for sitting furniture

Category: Under advice

Click [here](#) for more information

Warranty Information

- 2 Years

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
CHP227640FA	Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Soft Touch Fabric	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm
CHP227640LT	Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Leather	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm
CHP227640MS	Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Microsuede	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm
CHP227640VL	Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Healthcare Vinyl	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm
CHP227650FA	Raphael X CUSTOM Quattro Lift Recline Chair - SMALL - Soft Touch Fabric	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm
CHP227650LT	Raphael X CUSTOM Quattro Lift Recline Chair - SMALL - Leather	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm
CHP227650MS	Raphael X CUSTOM Quattro Lift Recline Chair - SMALL - Microsuede	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm
CHP227650VL	Raphael X CUSTOM Quattro Lift Recline Chair - SMALL - Healthcare Vinyl	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
CHP227660FA	Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM Soft Touch Fabric	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm
CHP227660LT	Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM - Leather	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm
CHP227660MS	Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM - Microsuede	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm
CHP227660VL	Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM - Healthcare Vinyl	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm
CHP227670FA	Raphael X CUSTOM Quattro Lift Recline Chair - LARGE Soft Touch Fabric	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm
CHP227670LT	Raphael X CUSTOM Quattro Lift Recline Chair - LARGE - Leather	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm
CHP227670MS	Raphael X CUSTOM Quattro Lift Recline Chair - LARGE - Microsuede	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm
CHP227670VL	Raphael X CUSTOM Quattro Lift Recline Chair - LARGE - Healthcare Vinyl	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm
CHP227680FA	Raphael X CUSTOM Quattro Lift Recline Chair - MAXI- Soft Touch Fabric	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm
CHP227680LT	Raphael X CUSTOM Quattro Lift Recline Chair - MAXI- Leather	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm
CHP227680MS	Raphael X CUSTOM Quattro Lift Recline Chair - MAXI- Microsuede	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
CHP227680VL	Raphael X CUSTOM Quattro Lift Recline Chair - MAXI - Healthcare Vinyl	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm
CHP229200	Classic Backrest - Petite - Golden Grain	-	-	-	-	-	-	-
CHP229201	Classic Backrest - Petite - Marina Grey	-	-	-	-	-	-	-
CHP229202	Classic Backrest - Petite - Nordic Stone	-	-	-	-	-	-	-
CHP229210	Classic Backrest - Small - Golden Grain	-	-	-	-	-	-	-
CHP229211	Classic Backrest - Small - Marina Grey	-	-	-	-	-	-	-
CHP229212	Classic Backrest - Small - Nordic Stone	-	-	-	-	-	-	-
CHP229220	Classic Backrest - Medium - Golden Grain	-	-	-	-	-	-	-
CHP229221	Classic Backrest - Medium - Marina Grey	-	-	-	-	-	-	-
CHP229222	Classic Backrest - Medium - Nordic Stone	-	-	-	-	-	-	-
CHP229230	Classic Backrest - Large - Golden Grain	-	-	-	-	-	-	-
CHP229231	Classic Backrest - Large - Marina Grey	-	-	-	-	-	-	-
CHP229232	Classic Backrest - Large - Nordic Stone	-	-	-	-	-	-	-
CHP229240	Classic Backrest - Maxi - Golden Grain	-	-	-	-	-	-	-
CHP229241	Classic Backrest - Maxi - Marina Grey	-	-	-	-	-	-	-
CHP229242	Classic Backrest - Maxi - Nordic Stone	-	-	-	-	-	-	-