

Compression Armsleeve - Grip Top



- Designed for security and support all day long
- Medically graduated compression for lymphedema symptoms
- High modulus inlay yarns to maintain compression without garment fatigue

 Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

| Code | Display Name | Elbow Circumference |
|-----------|---|---------------------|
| FBS310240 | Compression Armsleeve with Grip Top - Large L2 - Long - Beige | 260 - 320 mm |
| FBS310250 | Compression Armsleeve with Grip Top - Medium M2 - Long - Beige | 230 - 290 mm |
| FBS310260 | Compression Armsleeve with Grip Top - Medium M1 - Short - Beige | 230 - 290 mm |
| FBS310270 | Compression Armsleeve with Grip Top - Small S2 - Long - Beige | 200 - 260 mm |
| FBS310280 | Compression Armsleeve with Grip Top PLUS - Medium M4 - Long - Beige | 290 - 370 mm |