

## Kings Button Hook



- Assists one handed dressing
- Pass through the button hole and over the button to grip the thread
- Button pulled through with twisting motion
- Shaped stainless steel wire with plastic handle
- Wire Length: 97 mm
- Weight: 75 g

### Warranty Information

- 1 Year

- ✓ Contracted DVA product
- ✓ Contract approved

### ndis

#### Item Code:

05\_090900046\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

| Code      | Display Name      | RRP  | Box Quantity | Total Length |
|-----------|-------------------|------|--------------|--------------|
| DLD269300 | Kings Button Hook | 26.8 | 1            | 120 mm       |