

Kerry Hammock Sling - Mesh



- Attaches to hoist to enable a carer to lift and manoeuvre a less mobile user
- Provides excellent back and head support for patients requiring support during transfers
- Quick and easy to fit
- Ideal for non weight bearing patients
- Colour coded attachment and positioning straps
- Machine washable
- FIM Scores: 1 and 2



Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 1 Year

Code	Display Name	Safe Working Load	Material
LSS383200	Hammock Sling - Small	250 kg	Mesh
LSS383300	Hammock Sling - Medium	250 kg	Mesh
LSS383400	Hammock Sling - Large	250 kg	Mesh