

Mangar Eagle - Lift Cushion



The Eagle lifting cushion will both sit up and lift a fallen person. It has a small footprint and can be used in small spaces such as a bathroom. The small footprint design allows for the Eagle to be used in confined spaces, where traditional lifting equipment may not fit. Ideal for use within a home environment.



- Inflate swiftly using the Airflo Compressor and intuitive hand control—it is as easy as a touch. Minimal training is required, making it accessible for users of all levels. Whether indoors or outdoors, its versatility knows no bounds. Use it independently or with assistance.
- Easy to clean using a proprietary, non-abrasive mild liquid cleaner, simply deflate, fold, and stow away in the provided carry bag – ideal for infection control management.
- Designed with safety in mind. Capable of lifting individuals weighing up to 220kg from the floor, it provides unparalleled rigidity and support to the back.
- Perfect for supported living apartments and home care settings, it offers peace of mind with every lift.
- Experience a lift in both safety and dignity. By reducing the risk of musculoskeletal injuries for caregivers, it is a game-changer in care routines.
- Individuals with cognitive impairments find comfort in the Eagle's gentle lifting mechanism, while the added backrest support provides comfort and peace of mind for carers and users alike.

Code	Display Name	RRP	Box Quantity
------	--------------	-----	--------------

LSS403800	Mangar Eagle - Lift Cushion - 220kg	7019	1
-----------	-------------------------------------	------	---

Code	Display Name	RRP	Box Quantity
------	--------------	-----	--------------