

Walking Stick Tip



- The Hugo® Claw™ Standing Cane Tip lets your cane stand on its own!
- Lets your cane stand on its own! Stability in every direction.
- 57% lighter than the QuadPod, weighing only 0.2 lb (92 grams).
- 360° grip: great for folding canes!
- Tool-free installation compatible with all canes with 3/4" (1.9 cm) cane tubing.
- There are no reviews for this product.

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Walking, manipulated by one arm

Category: Prescribed

Click [here](#) for more information

Code

Display Name

WAS747429 Walking Stick Tip - Claw - Mini Stable Tip