

## Height Adjustable Toilet Seat Raiser



The Height Adjustable Toilet Seat Raiser helps you safely get on and off the toilet. Its sturdy arms provide stable support and the easy-to-clean frame sits snugly and subtly over the existing bowl.

- Height Adjustable to 50mm, 100mm and 150mm.
- Swing-away arms for convenient side access (43cm width between armrests).
- Sturdy plastic design with a Safe Weight Load of 120kgs.
- Easy to clean plastic surface.

✔ Contract approved

### Support @ Home

**Function:** Toileting

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Total Width	Total Height	Safe Working Load
BTT157950	Adjustable Toilet Seat Raiser w/ Arms	430 mm	50 - 150 mm	120 kg