

## INVISA-BEAM Free-Standing Monitor



- INVISA-BEAM is an early warning monitoring device that gives positive detection for falls prevention
- The MPM (Multi-Purpose Monitor) free standing monitor is placed at the head of the bed to monitor the side of the bed
- Watch [informational video here](#)
- Uses invisible beams to define the edge of the bed
- Precise pencil point beam is tuned to the length of a bed to monitor exit and entry from the bed only
- Can be easily transferred from bed to bed
- Hygienic wipe clean smooth surface for easy infection control
- Constantly monitors the bed patient/resident
- Can be utilised with all Invisa-Beam wireless alarms as well as connected to nurse call.
- Alerts nurse/carer immediately if patient/resident tries to leave their bed
- It is mains power operated
- The height is adjustable for lowered beds

✔ Contract approved

### Support @ Home

**Function:** Signalling, alarming and localisation

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 12 Months
- 24 Months Electronics

Code	Display Name
BEA007222	MPM Bed Monitor w/ INVISALERT

Code	Display Name
BEA007220	MPM Bed Monitor w/ Mono & Stereo Plug
BEA007221	MPM Bed Monitor w/ Remote Alarm