

Coway BA-08 Electric Bidet



The Coway BA08 Bidet is a versatile and feature-rich bidet designed to suit the needs of the entire family. It is an ideal solution for personal health and hygiene, providing a range of advanced features for comfort and effectiveness. With innovative technology and user-friendly design, the BA08 ensures a superior cleaning experience for everyone.

- **One-Touch Button:** Simplifies operation with a single button for wash and dry functions.
- **Adjustable Spray Options:** Choose between narrow, medium, and wide spray settings for a customised wash experience.
- **Air+ Button:** Provides a gentler yet more effective wash.
- **Dual Nozzle System:** Separate nozzles for rear and front washes reduce the risk of cross-contamination.
- **Extra-Large Remote Control:** Water-resistant, with an LCD display and braille buttons, it can be handheld or wall-mounted.
- **Optional One Touch Button Tamper Proof Holder:** Available for added convenience and security.
- **Ideal for Various Conditions:** Suitable for individuals with conditions like Dementia, Alzheimer's, Parkinson's Disease, Multiple Sclerosis, Limited Mobility, Arthritis, and more.
- **Custom Seat Colours:** Coloured seats available to assist with colour differentiation, especially for individuals with dementia or Alzheimer's. Contact us for more details.
- **Installation Requirements:** Ensure access to seat bolts from underneath the toilet. Power lead on the RHS and water pipe on the LHS—these cannot be modified.

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_093318082_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Warranty Information

- 2 Years Limited - 1 Year Standard Manufacturer's Warranty, 1 Year Complimentary Parts & Labour Warranty

Code	Display Name	Total Length	Total Width	Total Height	Unit Weight
BTT150110	Coway BA-08 Bidet - Warm Water - Elongated	530 mm	525 mm	142 mm	5.7 kg - 6.3 kg
BTT150130	Coway BA-08 Bidet - Warm Water - Non Elongated	500 mm	525 mm	142 mm	5.7 kg - 6.3 kg