

## ORTHOBRACE OBWRI1 ORTHOWRIST



- Extra-long forearm design
- 3 straps provide optimal support\*
- Volar, dorsal and lateral aluminium stays
- Soft inner lining Application:
- Buckle fractures
- Mild, stable distal radius and wrist fractures
- Support following cast removal
- Wrist sprains and strains

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name
FBS308270	ORTHOBRACE OBWRI1 ORTHOWRIST - LEFT - XS
FBS308280	ORTHOBRACE OBWRI1 ORTHOWRIST - RIGHT - XS
FBS308290	ORTHOBRACE OBWRI1 ORTHOWRIST - LEFT - SML
FBS308300	ORTHOBRACE OBWRI1 ORTHOWRIST - RIGHT - SML
FBS308310	ORTHOBRACE OBWRI1 ORTHOWRIST - LEFT - MED
FBS308320	ORTHOBRACE OBWRI1 ORTHOWRIST - RIGHT - MED
FBS308330	ORTHOBRACE OBWRI1 ORTHOWRIST - LEFT - LGE

**Code**

**Display Name**

FBS308340 ORTHOBRACE OBWRI1 ORTHOWRIST - RIGHT - LGE