

Aidacare Aspire Nova Standing Lifter



The Nova Standing Lifter combines a lightweight aluminium frame, high safe working load lift solution for professional care environments. Its modern design and aluminium frame make it ideal for facilities seeking a safer, more efficient way to perform standing lifts in tight-room access for variable user sizes. The Nova is competitively priced, with unique features to cater for the demands of a commercial healthcare environment.



- **Bonus Second TiMotion Battery & External Charging Cradle** supplied doubles the operational capacity and reduces lifter downtime
- **39.6kg Aluminium Frame with Compact Base Design** - Supports a 250kg safe working load, reduces carer effort in facilities with crowded space
- **Motorbike-Style Handles** - Enhances user confidence and stability during transfers
- **Colour-Coded Sling Locks** - Patented Design enables quick visual checks for safer, more confident transfers
- **Emergency Lowering & Stop Functions** - Built-in safety for power or system failures
- **Over 1000mm Lifting Range** - Supports a broad spectrum of user mobility needs, ensures safe, efficient clinical transfers
- **Durable Leg Shrouds** - Protects components when going under beds, maintains a professional look and are easily replaceable
- **Adjustable Knee Pads & Leg Strap** - Comfort-focused design with infection control surfaces
- **Smart Handset** - With battery levels and service status indicators, supports efficient product servicing
- **Powered Leg Opening & Low-Mounted Footplate** - Allows smooth positioning around wheelchairs and chairs, improving both access and safety
- **Interchangeable Clip Boom** - Optional accessory allows compatibility with

both loop and clip slings for flexible lift handling

- **100mm Twin Front Castors & 125mm Twin Rear Castors** - Smaller front castor enhances underbed clearance and manoeuvrability, with rear locks for added control
- **Broad Sling Compatibility** - Works with Aspire and most third-party slings

Warranty Information

- 5 Years - Frame
- 3 Years - Electronics
- 1 Year - Battery, Handset & Cables

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range
LSS390652	Aidacare Aspire Nova Aluminium Standing Lifter	4595	1	1100 mm	682 - 1138 mm	1095 - 1880 mm	250 kg	Aluminium	39.6 kg	692 - 1770 mm
LSS390660	Nova Standing Lifter Interchangeable Clip Boom	799	1	-	-	-	-	-	-	-