

## Rolyan - Breathoprene Wrist Splint

Experience superior wrist support with the Breathoprene Wrist Splint - the one you'll always want to wear!

✔ Contract approved

### Key Features:

- Breathable, stretchy 3.2mm Breathoprene material
- Removable metal stays and D-ring wrist strap for added support
- Latex free, easy to wash and maintain

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information



### Enhancements:

- Stronger dorsal and palmar stays; longer splint length
- Wider forearm sizing for a better fit
- Longer, wider hand and forearm straps for improved comfort

### Care Instructions:

- Remove metal stays before washing
- Hand wash in mild detergent, air dry only

Code	Display Name	Wrist Circumference
FBS301417	Rolyan - Breathoprene Wrist Splint - Sml - Left	14 - 16.5cm
FBS301418	Rolyan - Breathoprene Wrist Splint - Sml - Right	14 - 16.5cm
FBS301419	Rolyan - Breathoprene Wrist Splint - Med - Left	16.5 - 18.4cm

<b>Code</b>	<b>Display Name</b>	<b>Wrist Circumference</b>
FBS301420	Rolyan - Breathoprene Wrist Splint - Med - Right	16.5 - 18.4cm
FBS301421	Rolyan - Breathoprene Wrist Splint - Lge - Left	18.4 - 20.3cm
FBS301422	Rolyan - Breathoprene Wrist Splint - Lge - Right	18.4 - 20.3cm
FBS301423	Rolyan - Breathoprene Wrist Splint - XL - Left	20.3 - 22.2cm
FBS301424	Rolyan - Breathoprene Wrist Splint - XL - Right	20.3 - 22.2cm