

K Care Hand Made Commode Seat - Open



- Superior quality white marine grade vinyl provides durability and comfort
- Vinyl is wrapped, stapled and sealed, making cleaning easier with no seam on the inside of the aperture
- Comes with plastic seat runners for easy attachment and removal for cleaning
- Seat height is 55 mm

Warranty Information

- 1 Year



Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Material	Seat Depth	Seat Width
BTC073900	Open Front Commode Seat - Hand Made - 400 mm Wide	Vinyl and Foam	430 mm	400 mm
BTC074400	Open Front Commode Seat - Hand Made - 440 mm Wide	Vinyl and Foam	430 mm	440 mm
BTC062800	Open Front Commode Seat - Hand Made - 500 mm Wide	Vinyl and Foam	430 mm	500 mm
BTC062900	Open Front Commode Seat - Hand Made - 550 mm Wide	Vinyl and Foam	430 mm	550 mm
BTC063000	Open Front Commode Seat - Hand Made - 600 mm Wide	Vinyl and Foam	430 mm	600 mm
BTC063100	Open Front Commode Seat - Hand Made - 650 mm Wide	Vinyl and Foam	430 mm	650 mm
BTC074000	Left Side Open Commode Seat - Hand Made - 440 mm Wide	-	430 mm	440 mm

Code	Display Name	Material	Seat Depth	Seat Width
BTC074500	Right Side Open Commode Seat - Hand Made - 440 mm Wide -		430 mm	440 mm