

K Care Hand Made Commode Seat - Open



- Superior quality white marine grade vinyl provides durability and comfort
- Vinyl is wrapped, stapled and sealed, making cleaning easier with no seam on the inside of the apperture
- Comes with plastic seat runners for easy attachment and removal for cleaning
- Seat height is 55 mm

Warranty Information

• 1 Year

Code	Display Name	Seat Depth	Seat Width	Material
BTC073900	Open Front Commode Seat - Hand Made - 400 mm Wide	430 mm	400 mm	Vinyl and Foam
BTC074400	Open Front Commode Seat - Hand Made - 440 mm Wide	430 mm	440 mm	Vinyl and Foam
BTC062800	Open Front Commode Seat - Hand Made - 500 mm Wide	430 mm	500 mm	Vinyl and Foam
BTC062900	Open Front Commode Seat - Hand Made - 550 mm Wide	430 mm	550 mm	Vinyl and Foam
BTC063000	Open Front Commode Seat - Hand Made - 600 mm Wide	430 mm	600 mm	Vinyl and Foam
BTC063100	Open Front Commode Seat - Hand Made - 650 mm Wide	430 mm	650 mm	Vinyl and Foam
BTC074000	Left Side Open Commode Seat - Hand Made - 440 mm Wide	430 mm	440 mm	-



Compliance

AS/NZS 3973:2009



Code	Display Name	Seat Depth Seat Width	Material	
------	--------------	-----------------------	----------	--

BTC074500 Right Side Open Commode Seat - Hand Made - 440 mm Wide 430 mm 440 mm -