

K Care Hand Made Commode Seat - Open



- Superior quality white marine grade vinyl provides durability and comfort
- Vinyl is wrapped, stapled and sealed, making cleaning easier with no seam on the inside of the aperture
- Comes with plastic seat runners for easy attachment and removal for cleaning
- Seat height is 55 mm



Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year

Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width	Material
BTC073900	Open Front Commode Seat - Hand Made - 400 mm Wide	430 mm	400 mm	Vinyl and Foam
BTC074400	Open Front Commode Seat - Hand Made - 440 mm Wide	430 mm	440 mm	Vinyl and Foam
BTC062800	Open Front Commode Seat - Hand Made - 500 mm Wide	430 mm	500 mm	Vinyl and Foam
BTC062900	Open Front Commode Seat - Hand Made - 550 mm Wide	430 mm	550 mm	Vinyl and Foam
BTC063000	Open Front Commode Seat - Hand Made - 600 mm Wide	430 mm	600 mm	Vinyl and Foam
BTC063100	Open Front Commode Seat - Hand Made - 650 mm Wide	430 mm	650 mm	Vinyl and Foam
BTC074000	Left Side Open Commode Seat - Hand Made - 440 mm Wide	430 mm	440 mm	-

Code	Display Name	Seat Depth	Seat Width	Material
BTC074500	Right Side Open Commode Seat - Hand Made - 440 mm Wide	430 mm	440 mm	-