

## NuboMed X3000 Lite Mobile Nursing Cart



NuboMed

- **Material:** Constructed from a combination of ABS and aluminium alloy
- **Integrated Design:** The handle and keyboard tray are seamlessly integrated from a single piece of anti-microbial ABS material, ensuring a gap-free structure inclusive of a spill-proof groove design
- **Expanded Work Surface:** Pulls out from the rear handle for convenient use
- **Double-Pillar Design:** Cabinet securely mounted to the dual aluminium alloy pillars
- **Monitor Mount:** Height adjustable from 0-150mm, with 90° left and right swivel, 30° upward tilt, and portrait-to-landscape rotation
- **Button:** PC power on/off button and a screen backlight on/off switch
- **Drawer Housing:** ABS drawer housing with lock, equipped with four configurable aluminium alloy drawers (1U+2U+1U+2U). 5kg SWL per drawer and dividers are included
- **Base:** Designed to prevent tripping and made from a single-piece of aluminium alloy with an antibacterial powder coating finish
- **Castors:** Equipped with 2x 4" lockable and 2x 4" tracking castors
- **Ideal For:** Daily Nursing/Data Recording



### Warranty Information

- 2 Years

Code	Display Name
WAT908180	NuboMed - X3000 - Lite Mobile Nursing Cart

