

2 Section Bobath Examination Couch



- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Minimum height 50 cm to maximum height 90 cm in 30 seconds
- Adjustable backrest featuring double stainless steel struts, angle adjustable from 0 to +70°
- 3 breath holes with plugs in head section
- Adjustable foot for uneven surfaces
- 7.5 cm castors for safety and mobility
- Retractable locking castors



ndis

Item Code:

05_093312081_0103_1_2

Walking Aids -

Composite:05_415_0103_1_2

Click [here](#) for more information

Warranty Information

- 5 Years Frame
- 3 Years Electrics

Compliance

IEC 60601-2-52:2009

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Colour
WAT821570	Exam Couch - 2 Section Bobath Neurological Hi/Lo - Navy Blue	1890 mm	1200 mm	500 - 900 mm	250 kg	Steel/Vinyl/Foam	Navy