

DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB



- Dedicated for either right or left hand use
- Ergonomic, flexible lightweight design
- Anatomically designed, provides tailored individual fit
- Easy and simple application via adjustable straps and 'Stop Strap' system
- Incorporates two semi-rigid dorsal stays and an aluminum palmar stay for overall wrist immobilisation
- Advanced Ozone™ material, promotes air flow
- Breathable in general use, providing good skin care and around-the-clock comfort in different heat and humidity conditions Application:
- Wrist sprains and trauma
- Tendinopathy
- Pre and post-operative immobilization

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS305950	DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB - LEFT - XS
FBS305960	DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB - LEFT - SML
FBS305970	DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB - LEFT - MED
FBS305980	DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB - LEFT - LGE
FBS305990	DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB - LEFT - XL
FBS306000	DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB - RIGHT - XS
FBS306010	DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB - RIGHT - SML

Code	Display Name
FBS306020	DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB - RIGHT - MED
FBS306030	DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB - RIGHT - LGE
FBS306040	DONJOY 82-9617/8 RESPIFORM WRIST AND THUMB - RIGHT - XL