

Sigvaris - Traditional Sleeves



Arm compression sleeves are essential for the prevention, treatment, and care of lymphedema in the upper extremities, playing a vital role in maintenance therapy.

- **Skin-friendly:** Featuring a large-mesh fabric structure, these sleeves offer high breathability, ensuring comfort during wear.
- **Pressure stability:** Constructed with natural rubber, these sleeves provide a high recovery force, ensuring consistent pressure for effective compression therapy.
- **Wear-resistant:** Proven to be up to 6 times better in abrasion resistance compared to jeans, these sleeves offer long-lasting durability, making them suitable for daily wear.
- Traditional Sleeve w/ Mitten provide extra coverage and play a vital role in maintenance therapy.

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name
PTS642415	Sigvaris - Traditional Sleeves (specify size/colour)
PTS642420	Sigvaris - Traditional Sleeve w/ Mitten (specify size/colour)