

## Standing Lifter/Hoist 200kg - Sling separate



The Aidacare Aspire 200S is a robust and durable steel standing lifter capable in the most demanding facility environments



- Ergonomic Handles. Provide multiple grip points to suit carers of different heights and strengths
- Ergonomic Handset. The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Powered Leg Function. Opens and closes the legs for lifting stability or accessibility
- Sturdy steel frame. with industry leading Linak electronics and warranty provision
- Safety. Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- Adjustable Padded Legrests. Provide comfort and support during transfers & can be adjusted to suit the patient
- Intelligent Diagnostics. State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance

\*\*\* Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

**Code**

**Display Name**

H\_LSS390650 Standing Lifter/Hoist 200kg - Sling separate

