

Bedside Cabinet - 2 Drawers + 1 Dump Drawer + Lock



- This Bedside Cabinet features two standard drawers and one larger dump drawer at the bottom
- Lockable top drawer
- Easy to clean finish
- Features top edging to prevent items from falling off
- Colour options available



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height	Colour
BEA828113	Bedside Cabinet - 2 Drawers + 1 Dump Drawer + Lock (Beech)	460 mm	450 mm	810 mm	Beech
BEA828114	Bedside Cabinet - 2 Drawers + 1 Dump Drawer + Lock (White)	460 mm	450 mm	810 mm	White
BEA828115	Bedside Cabinet - 2 Drawer + 1 Dump Drawer + Lock (Sublime Teak)	460 mm	450 mm	810 mm	Sublime Teak
BEA828116	Bedside Cabinet - 2 Drawer + 1 Dump Drawer + Lock (Milano Walnut)	460 mm	450 mm	810 mm	Milano Walnut