

## Malvestio Paediatric Ward Bed



- Three function electrically adjustable hospital bed designed with growing people in mind
- Colourful side rails and decals to meet interior design themes and add to a more cheerful patient environment
- Capacity to lengthen at the press of a button to transition from child to young adult size proportions
- Biocote antibacterial surface treatment for improved infection control
- Durable blow molded side rails with gas spring dampeners for smooth operation
- Onboard battery provides ward mobility and emergency backup
- Integrated control panel in foot end section for carer and parent access
- Backrest angle 0 to 70 degrees
- Knee bend angle 0 to 36 degrees
- Trendelenburg -15 to 15 degrees

✓ Contract approved

### Compliance

IEC 60601-2-52:2009

## Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
PED490110	Mattress - Malvestio Paediatric Ward Bed	1580 mm	790 mm	140 mm	150 kg
PED490100	Malvestio Paediatric Ward Bed	1830 - 2180 mm	890 mm	550 - 950 mm	150 kg