

## Bed/Chair Raiser



- Raisers are placed on the floor under the legs of beds or lounge chairs to raise their height
- Ideal for those recovering from hip replacements, knee operations or as a result of ageing

### Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Height adjustment of furniture

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_181212202\_0103\_1\_2

**Minor Bed – Composite:**

05\_428\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Total Height
BEA009700	Bed / Chair Raiser - Blocks - Wooden 75mm (Pkt 4)	75 mm

Code	Display Name	Total Height
BEA009600	Bed/Chair Raiser Cones - Plastics - Pack of 4 - 140mm High	140 mm