

DermaSaver Shin Tube



- Two layers of soft, lightweight MicroSpring Textile form a loose-fitting tube
- Designed to protect fragile skin from pressure, friction and minor traumas that result in tears, bruising and abrasions
- Sizing: Simply measure the circumference at the widest part of the calf. For borderline measurements, choose the larger size for X comfort

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Calf/Arm Circumference
PTA533322	DermaSaver - Shin Tube - X Sml	18-23cm
PTA533322S	DermaSaver - Shin Tube - Sml	23-30cm
PTA533322M	DermaSaver - Shin Tube - Med	30-38cm
PTA533322L	DermaSaver - Shin Tube - Lge	38-48cm
PTA533322XL	DermaSaver - Shin Tube - XL	48-58cm