

Freedom Bedside Commode - JBA500-C



- Allows for toileting in the bedroom for users unable to travel between the bedroom and bathroom
- Removable padded backrest
- Padded armrest
- Hinged flip-up padded seat with spill proof upholstery
- Anti-slip suction rubber tips
- Height adjustable
- Removable durable plastic seat
- Includes pail and lid
- Available in beige

✓ Contracted DVA product

✓ Contract approved

Compliance

AS/NZS 3973:2009

Warranty Information

- 1 Year

freedom
HEALTHCARE

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Safe Working Load	Unit Weight	Seat Depth	Seat To Floor Height
BEC046405	Freedom Bedside Commode - With Luxury Padded Upholstery	427	1	600 mm	620 mm	160 kg	10.1 kg	490 mm	500 - 650 mm