

## ASPIRE REHAB RX JUNIOR WHEELCHAIR



✓ [Contract approved](#)

i **ndis**

Item Code: 05\_122203115\_0106\_1\_1

Wheelchair Manual – Composite:

05\_417\_0105\_1\_2

Click [here](#) for more information

### Accessories

Click on the cart icon to display the range of optional accessories



- **The Rehab RX Junior** - is designed for the increased postural support and pressure care of children, via easy repositioning and advanced adjustability.
- **Tilt (26°) and backrest recline (30°)** – this facilitates weight re-distribution and the preservation of skin integrity. It also optimises user comfort, decreases user fatigue and assists with head control.
- **Highly adjustable** - depth and width adjustable seat sizing. Height adjustable and removable armrests. Adjustable headrest.
- **Contoured, built-up backrest and seat cushion** - the backrest has a gentle lumbar build up and moderate lateral supports; together with the abductor and adductor build-up at the seat cushion, midline positioning is optimised.
- **Breathability** - support surfaces offer multi-layered foam with breathable multi-stretch covers - for pressure re-distribution, comfort and to enhance skin microclimate.
- **Advanced articulating legrests** – closely mimic the normal movement at the knee upon elevation, allowing the thigh to maintain contact with the seating surface, ensuring maximum support and skin contact with the cushion.
- **Crash Tested** – to AS/NZ3616.196

## Warranty Information

- 2 Years Frame
- 1 Year Components

Aspire Rehab RX Junior - 320-370mm Wide  
MWS449630  
EACH

Enquire EACH

SPECIFICATIONS						
Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load
MWS449630	Aspire Rehab RX Junior - 320-370mm Wide	320 - 370 mm	320 - 370 mm	460 - 710 mm	400 - 475 mm	180 kg

1300 216 898