

## iCare Contour ActiveX Pillow



- The Care ActiveX™ Contour neck pillow is made up of a high-density ActiveX™ core in a contoured pillow shape, a cooling gel insert, and a soft puff quilt style cotton cover. This Contour ActiveX™ Pillow is recommended for back and neck pain sufferers
- Contoured shape to support neck & shoulders
- Breathable vant fabric to increase airflow
- Cooling gel insert
- Best suited to back & side sleepers
- Suitable for pain management and alignment
- Contour design increasing pressure redistribution to head and neck
- Gel insert assists with temperature regulation

✔ Contract approved

### Support @ Home

**Function:** Body stabilisation

**Category:** Under advice

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height
BEP051060	iCare Contour ActiveX Pillow	600 mm	400 mm	90 - 110 mm