

Deluxe Paediatric Cot



- Paediatric Cot Deluxe features a stainless steel or epoxy coated frame
- Mechanisms to operate drop sides and ends located outside of the cot at mattress base for added safety
- Standard foam 50 mm mattress
- Four castors - two lockable
- Can be tilted at Trendelenburg angle, up to 12 degrees above the horizontal
- Meets AS 2130



Compliance

No applicable compliance standard

Warranty Information

- 3 Years

Code	Display Name	Total Length	Total Width
WAT831300	Paediatric Cot Deluxe - Stainless Steel	1440 mm	760 mm
WAT831800	Paediatric Cot Deluxe - White	1440 mm	760 mm
WAT830000	Paediatric Cot Chart Holder	-	-