

## Phoenix Chair - Paediatric



- The Paed Phoenix is the most supportive and functional of the Seating Matters chairs for children
- The Paed Phoenix chair provides pressure distribution due to the multi-stretch material and reflex foam in the seat cushion and offers a great deal of support with the tilt in space mechanism and adjustable supports
- The head support accommodates those with kyphotic postures, hyperextension, or those who require upper body and head support
- The angle of the back can be adjusted to accommodate the patient and the elevating legrest can be raised and lowered, also increasing patient comfort
- The tilt in space will allow all angles on the chair to remain the same while tilting the patient back, allowing their weight to be spread throughout the body and reduce pressure on a given point
- The push handle and castors allow the patient to be transported in the chair. It is fully hoist accessible which reduces moving and handling significantly
- The removable arms facilitate width adjustment and also side transfers
- The arm height, seat depth and width as well as the footplate height and angle are all fully adjustable
- Effective for children with: Risk of pressure ulcers, Reduced upper body control, Kyphosis, Scoliosis, Muscle contractures & Swelling

### Warranty Information

- 7 years frame
- 1 year upholstery

 Contract approved

### Support @ Home

**Function:** Sitting furniture

**Category:** Under advice

Click [here](#) for more information

### NDIS

**Item Code:**

05\_180939182\_0105\_1\_2

Click [here](#) for more information

### Compliance

ISO13934-1 / ISO13937-1 / BS6807 / BS7175

### Accessories

Click on the cart icon to display the range of optional accessories

Code	Display Name	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Backrest Height
CHP197800	Chair Pressure Relief - Phoenix Chair - Paediatric - Manual - 160kg	1220 mm	160 kg	57 kg	290 - 380 mm	250 - 400 mm	540 mm	680 - 830 mm