

## Aidacare Aspire 200 Steel Patient Lifter



A robust workhorse, the Aidacare Aspire 200 is designed to meet the rigorous demands of professional and facility environments

- **Ergonomic Handles** Provide multiple grip points to suit carers of different heights and strengths
- **Ergonomic Handset** The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- **Powered Leg Function** Opens and closes the legs for lifting stability or accessibility
- **Quick Release Yoke Attachment** One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety
- **Sturdy steel frame** with industry leading Linak electronics and warranty provision
- **Safety** Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- **Intelligent Diagnostics** State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- **Tente Casters:** 2 x 75 mm front castors for maximum under bed clearance and accessibility. 2 x 100 mm rear locking castors for safety and manoeuvrability.
- Item Weight: 43.9 kg
- Tested and compliant to AS/NZS ISO 10535:2011 (the Australian standard for hoists)

### Warranty Information

- Frame: 5 Years

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Lifting persons

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_123604167\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

### Accessories

Click on the cart icon to display the range of optional accessories

- Electronics: 3 Years
- Battery and Handset: 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range	Electric Leg
LSS390620	Aidacare Aspire 200 Patient Patient Lifter	1280 mm	665 - 980 mm	1940 mm	200 kg	Steel	43.9 kg	535 - 1725 mm	580 - 900 mm	True