

## Sleepform



- Sleepform is a vacuum mattress moulded to support the user as part of a 24 hour positioning program
- Sleepform System is designed for children, adolescents and adults with special needs
- Promotes symmetry, reduces pain, slows the development of deformities
- Designed for use at home or in hospital and respite care centres



**ndis**

**Item Code:**

05\_181224204\_0103\_1\_2

**Minor Bed – Composite:**

05\_428\_0103\_1\_2

Click [here](#) for more information

## Warranty Information

- 2 Years

**LECKEY®**

Code	Display Name
PED477010	Sleepform - Age 0-1Y - Cot - Size 1
PED477020	Sleepform - Age 1-5Y - Cot - Size 2
PED477030	Sleepform - Age 1-5Y - Bed - Size 3
PED477040	Sleepform - Age 5-18Y - Bed - Size 4