


## Lifecomfort Mattress - Pressure Reducing (1940 x 900 x 150) - Med Pressure Care



- Pressure reducing mattress
- Four layer foam support design for maximum comfort and support
- Full size memory foam overlay
- Foam side walls for extra stability on the edge of the mattress
- Two way waterproof stretch cover with welded seams
- PVC bottom cover
- Three sided zip for easy cover removal
- Welded zipper cover with waterproof flap for the best hygiene
- PUPPS 2 (2004) compliant



\*\*\* Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

Code	Display Name
H_BEM046610	Lifecomfort Mattress - Pressure Reducing (1940 x 900 x 150) - Med Pressure Care