

Insulin Pump Belt



The Glucology pump belt is fitted with a wide opening, allowing you to place and remove items easily without any hassle. Its mesh panel is designed for easy viewing and operation of your pump without needing to remove the belt or take it out of the pouch. The Glucology pump belt is made to last, using durable and resilient materials that will ensure many years of functionality. This stretchable waist pouch does not have any clips, zips, or attachments that can hurt or irritate your skin while in use.

- Carry your pump discreetly with the pump belt having a slim and compact profile to wear under your shirt, sweater or coat without looking bulky.
- Easy access to your pump and cable with the mesh pouch allowing you to view your pump or take it out of the pouch without having to take the belt off
- Made with high-quality stretch fabric that conforms to your waist with a discrete hook-and-loop dot providing extra protection for your pump
- Store small personal items in any of the pocket compartments on the belt
- Ideal for children and adults with 4 different sizes available

✓ Contracted DVA product

