

Shoalhaven Manual Recliner



- Air cell cushioning combined with slow release memory foam, for excellent pressure relieving capabilities
- Adjustable and removable head rest wings
- Compact postural design with tilt in space recline action and high leg rest elevation for reduced fluid retention
- Resilient water proof stretch fabric with anti-bacterial and non-allergenic properties
- Designed for use with standing machines and transfer hoists
- Low seat height for fall out safety and easier stand assist transfers
- Hospital quality 100 mm bearing glide locking casters

✓ Contracted DVA product

✓ Contract approved

Compliance

ISO13934-1 / ISO13937-1 / BS6807
/ BS7175

Warranty Information

- 1 Year

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Safe Working Load
CHP198950	Shoalhaven Manual Recliner	500 mm	520 mm	750 mm	130 kg