

DONJOY 82-0012 STABILAX HINGED ELASTIC KNEE BRACE



- Multidirectional Elasticity - Enhances proprioceptive stimulation through sensory response
- Anatomical 3-D Knit - Conforms to the leg and provides optimal fit
- Comfort in Flexion - Looser popliteal knit for unrestricted flexion
- Support after mild to moderate medial and lateral ligament injury
- Proprioceptive & neuromuscular stimulation through compression
- Return to sports activities

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS305740	DONJOY 82-0012 STABILAX HINGED ELASTIC KNEE BRACE - XS
FBS305750	DONJOY 82-0012 STABILAX HINGED ELASTIC KNEE BRACE - SML
FBS305760	DONJOY 82-0012 STABILAX HINGED ELASTIC KNEE BRACE - MED
FBS305770	DONJOY 82-0012 STABILAX HINGED ELASTIC KNEE BRACE - LGE
FBS305780	DONJOY 82-0012 STABILAX HINGED ELASTIC KNEE BRACE - XL
FBS305790	DONJOY 82-0012 STABILAX HINGED ELASTIC KNEE BRACE - XXL