

Aspire Comfy Sleeper Chair-Bed



Seamlessly transforming from a snug lounge chair into a supportive sleeping surface, this chair caters to the needs of loved ones who wish to stay close to their hospitalised family members. The chair is upholstered in medical-grade silicone vinyl, meeting rigorous infection control standards while delivering unparalleled comfort and functionality.

- **Infection Control:** Designed with antimicrobial foam and medical-grade silicone vinyl. It not only meets stringent cleaning requirements for acute care but also enhances the user's sleeping experience.
- **Versatile Transformation:** With an intuitive single-hand activation, the chair effortlessly transforms into a comfortable sleeping surface. This dual-purpose design conserves space and offers loved ones the opportunity to stay close to the patient in a hospital setting.
- **Effortless Assembly:** The chair can be quickly assembled and disassembled without the need for any tools, making it user-friendly and convenient.
- **Easy Mobility:** Equipped with recessed locking castors, the chair ensures easy positioning and accessibility, enhancing user experience.
- **Optimal Comfort:** The chair provides a completely flat sleeping surface, designed to maximize rest and relaxation, ensuring a comfortable stay for the user.



Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

NDIS

Item Code:

05_180915811_0103_1_2

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height	Material	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height
CHP198808	Aspire Comfy Sleeper Chair-Bed - Beige	2155 mm	970 mm	940 mm	Beige Vinyl	62.2 kg	Beige	610 mm	740 mm	610 mm

Code	Display Name	Total Length	Total Width	Total Height	Material	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height
CHP198809	Aspire Comfy Sleeper Chair-Bed - Navy	2155 mm	970 mm	940mm	Navy Vinyl	62.2 kg	Navy	610 mm	740 mm	610 mm