

Aspire Comfy Sleeper Chair-Bed





Seamlessly transforming from a snug lounge chair into a supportive sleeping surface, this chair caters to the needs of loved ones who wish to stay close to their hospitalised family members. The chair is upholstered in medical-grade silicone vinyl, meeting rigorous infection control standards while delivering unparalleled comfort and functionality.

- **Infection Control:** Designed with antimicrobial foam and medical-grade silicone vinyl. It not only meets stringent cleaning requirements for acute care but also enhances the user's sleeping experience.
- Versatile Transformation: With an intuitive single-hand activation, the chair effortlessly transforms into a comfortable sleeping surface. This dual-purpose design conserves space and offers loved ones the opportunity to stay close to the patient in a hospital setting.
- **Effortless Assembly:** The chair can be quickly assembled and disassembled without the need for any tools, making it user-friendly and convenient.
- **Easy Mobility:** Equipped with recessed locking castors, the chair ensures easy positioning and accessibility, enhancing user experience.
- Optimal Comfort: The chair provides a completely flat sleeping surface, designed to maximize rest and relaxation, ensuring a comfortable stay for the user.

Code	Display Name	•	Total Width	Total Height	Material	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height
CHP198808	Aspire Comfy Sleeper Chair-Bed - Beige	2155 mm	970 mm	940 mm	Beige Vinyl	62.2 kg	Beige	610 mm	740 mm	610 mm





Code	Display Name	•	Total Width	Total Height	Material	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height
CHP198809	Aspire Comfy Sleeper Chair-Bed - Navy	2155 mm	970 mm	940mm	Navy Vinyl	62.2 kg	Navy	610 mm	740 mm	610 mm