

DermaSaver Arm Tube Protection



- Soft and loose with ample airflow reaching the skin allowing for greater healing ability whilst offering protection
- Cushioned layers of MicroSpring Textile protect fragile skin from friction and pressure that can cause skin tears, bruising and abrasions
- Sizing: simply measure the circumference around the widest part of the forearm. For borderline measurements, choose the larger size to ensure comfort

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Calf/Arm Circumference
PTA533010	DermaSaver Arm Tube with Double Elbow - X Sml	13-18cm
PTA533020	DermaSaver Arm Tube with Double Elbow - Sml	18-25cm
PTA533030	DermaSaver Arm Tube with Double Elbow - Med	25-30cm
PTA533040	DermaSaver Arm Tube with Double Elbow - L	30-38cm
PTA533050	DermaSaver Arm Tube with Double Elbow - XL	38-46cm
PTA533308	DermaSaver - Arm Tube - X Sml	13-18cm
PTA533308S	DermaSaver - Arm Tube - Sml	18-25cm

Code	Display Name	Calf/Arm Circumference
PTA533308M	DermaSaver - Arm Tube - Med	25-30cm
PTA533308L	DermaSaver - Arm Tube - Lge	30-38cm
PTA533308XL	DermaSaver - Arm Tube - XL	38-46cm