

Aidacare Bariscale



- Australian made portable bariatric scale with a 400 kg capacity
- Stainless steel handrails for easy cleaning and sanitising
- Powder coated weigh pan for durability
- Open front design and robust handrails allow for the largest patients or residents to safely access the platform
- Base plate has been designed to minimise tripping hazards and increase accessibility with its super low clearance (38 mm) and tapered sides
- Advanced digital display with battery and power pack
- Nylon castors for easy transportation
- 4 x high quality industrial load cells provide precision weight readings up to 400 kg



Compliance

AS/NZS ISO 10535:2011



Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
LSS362005	Aidacare Bariscale	660 mm	660 mm	900 mm	400 kg