

## Diffuser Cushion

- Premium density pressure-relief cushion made from memory foam
- Ideal for wheelchairs or conventional seating



### Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

**ndis**

**Item Code:**

05\_181006190\_0105\_1\_2

**Pressure Care Cushion –**

**Composite:05\_429\_0103\_1\_2**

Click [here](#) for more information

### Compliance

ISO 7176-7:1996

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Pressure Care Rating	Material
GSS340100	Diffuser Cushion	160	1	-	-	-	-
GSS340103	Diffuser Cushion with Waterproof Cover	168	1	420 mm	440 mm	Low Risk	Memory Foam
GSS340104	Diffuser Cushion - Large	182	1	460 mm	490 mm	Low Risk	Memory Foam
GSS340105	Diffuser Cushion with Waterproof Cover - Large	182	1	460 mm	490 mm	Low Risk	Memory Foam
GSS340101	Waterproof Cover - Standard	83	1	420 mm	440 mm	Low Risk	-
GSS340106	Waterproof Cover - Large Size	89.5	1	460 mm	490 mm	Low Risk	-

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Pressure Care Rating	Material
------	--------------	-----	--------------	--------------	-------------	----------------------	----------